#### Estimado cliente

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Stand der Information: 09/10

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## Dear Customer

Thank you for having chosen a quality product from Dentaurum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at [www.dentaurum.de](http://www.dentaurum.de).

### 1. Manufacturer

Dentaurum GmbH & Co. KG  
Turnstraße 31 · 75228 Ispringen · Germany

### 2. General product description

The Hyrax® screw (Hygienic Rapid Expansion Screw) is a special screw for rapid palatal expansion. The retention legs of the palatal expansion screw are welded or soldered onto preformed metal bands or cap splints. Alternatively, they can be polymerized into acrylic splints. The bands or splints are cemented in the usual way.

### 3. Area of application

Used for rapid palatal expansion with a fixed appliance.

### 4. Fabricating a Hyrax® appliance using the example of preformed bands

Preformed bands are fitted on the maxillary left and right 1st bicusps and left and right 1st molars. These bands are transferred to a dental stone model using an impression or fitted to a dental stone model. The legs of the Hyrax® screw are then bent to fit the bands so that they can be soldered or welded into place (Fig. 1).

To prevent the legs from breaking during bending, we recommend using our universal bending tool (REF 611-122-00, Fig. 4). If one is not available, hold the legs against the body of the expansion screw, using strong flat pliers (Fig. 5).

The Hyrax® screw should „float“ freely at least 1-2 mm above the palate roof/gingiva (Fig. 2). Its posterior side should form a straight line with the mesio palatal cusps of the 1st molars. To avoid any leverage effect, it should be ensured that the occlusal level is parallel.

### 5. Information on use of Hyrax® screws

To increase strength, the bands can be joined on the buccal and/or lingual side with a piece of 1.0 mm remaloy® wire (REF 528-100-00). Laser weld or solder these wires above the legs. The palatal strengthening wire should project slightly beyond the bands and be supported by the adjacent teeth to allow better distribution of expansion forces.

Easy, problem-free bending of the retention legs of Hyrax® screws with the universal bending tool (REF 611-122-00), without the use of additional pliers (leverage principle).

To prevent the legs from breaking during bending, they should be held against the body of the expansion screw using strong flat pliers e.g. Maxi wire-bending pliers (REF 013-522-00).

When using the soldering technique please protect the body of the screw against excessive heat, e.g. with heat screen (REF 380-070-00).

### 6. Note for doctor and patient

Before the Hyrax® appliance is cemented into position, turn the spindle (screw) slightly 2-3 forth of a turn in the direction of the arrow to test the movability. Then turn back to the original position.

The finished Hyrax® appliance is cemented into position in the usual manner.

Following this, the screw is turned four times by 90° using the key supplied for this purpose. One full 360° turn of the spindle = 0.8 mm expansion.

The patient should turn the screw two times 90° both in the morning and evening (according to your instructions) until the desired turn is achieved.

**Important:** use only the blue safety key supplied, as this cannot accidentally be swallowed.

The safety key is designed to hold a maximum force of 35 N ± 3.5 kg. Improper use may lead to overload and consequently to failure thus potentially excluding any guarantee and liability.

Regular visits are recommended to monitor progress.

The Hyrax® appliance is secured against unwanted turning-back, after the desired degree of expansion is achieved, by fitting a ligature or blocking with acrylic. For stabilization the appliance should remain in the mouth for approximately 2-3 months.

A retainer should be worn following removal of the Hyrax® appliance.

### 7. Composition

The composition is included in the materials list, refer to the catalogue or [www.dentaurum.de](http://www.dentaurum.de).

### 8. Safety information

Always use the blue safety key supplied to prevent swallowing!

The product should not be used, if there is a known allergic reaction to one or more of the material components.

### 9. Information on single-use products

The product is intended for single use only. Reconditioning of a product that has already been used (recycling) and its reuse on a patient is not permitted.

### 10. Availability

Description	straight legs	bent legs	Quantity
Hyrax® Mini - 7 (Expansion 7 mm)	<b>602-800-10</b> <b>602-800-30</b>	<b>602-801-10</b> <b>602-801-30</b>	1 piece 10 pieces
Hyrax® Medium - 10 (Expansion 10 mm)	<b>602-808-10</b> <b>602-808-30</b>		1 piece 10 pieces
Hyrax® Maxi - 12 (Expansion 12 mm)	<b>602-811-10</b> <b>602-811-30</b>	<b>602-812-10</b> <b>602-812-30</b>	1 piece 10 pieces

with safety key

Safety key	<b>611-120-01</b>	10 pieces
Universal bending tool	<b>611-122-00</b>	1 piece

### 11. With reference to our quality standards

Dentaurum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

### 12. Explanation of symbols used on the label

REF	Order number
12345	Batch description
0483	CE Mark MP Class IIa
→ <a href="http://www.dentaurum.de">www.dentaurum.de</a>	Instructions for use available in the internet
	Not for reuse
	Codablock barcode
	Matrix code
	Green dot
Date of information: 09/10	
Subject to modifications	

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